



# SPECIAL, Very Limited Time Offer!

One Device yields 100% ESI Sample Intro. Efficiency & makes Great MALDI Crystals!



nanoLiter Lyte, \$1.9 k!

ASMS Booth 57 !



Improves sensitivity, 10, 20x !

Pipette /MS ,  
Handheld or robotic.




Our NEW US Army's APG-X1 Dispenser.





# Programmable, nanoLiter Wave.

Optional Cell Phone Control!

Fastest (msec), Most Efficient (100%), ESI/MALDI,  
Sample Introduction/Prep Device In The World, LITERALLY!





Positive ion profiles generated by individual drops. a. Total ion current. b.  $[M(NO_3)_2(CMPO)]^+$ . c.  $[M(NO_3)(CMPO)_2]^{2+}$ . Blue =  $Ce^{3+}$ , Green =  $Tb^{3+}$ , Red =  $Lu^{3+}$ .



## Same Device Makes Excellent MALDI, SIMS, LDI Crystals!

Increase Sensitivity for MALDI, SIMS & LDI by 10-100x Literally! Ask for peer review references!



← Lo 3D Printing ?!

Contact Drew Sauter,  
nanoLiter LLC, nanoLiter.com, 702-882-5413